
# Estrogen Receptor (ER-12)

 CELL MARQUE DATA SHEET

## Estrogen Receptor (ER-12)

English .....	2	Italiano .....	18
Česky.....	4	Norsk.....	20
Dansk .....	6	Polski.....	22
Nederlands taal .....	8	Português .....	24
Français .....	10	Slovenski .....	26
Deutsch .....	12	Español .....	28
Ελληνικές .....	14	Svensk .....	30
Magyar .....	16		

### References:

1. Dabbs DJ et al. Diagnostic Immunohistochemistry 2002 Churchill Livingstone
2. Kell DL et al. Applied Immunohistochemistry 1(4): 275-81, 1993
3. Leong ASY et al. Applied Immunohistochemistry 1(4): 282-288, 1993
4. Tesch M et al. Am J Clin Pathol 1993; 99:8-12
5. Francesco AM et al. Applied Immunohistochemistry 2 (3): 157-163, 1994
6. Silvio M et al. Applied Immunohistochemistry 3(2): 85-90, 1995
7. Bejar J et al. Arch Pathol Lab Med 1998 Apr;122(4):346-52

**Presentation:**

Anti-Estrogen Receptor is a mouse monoclonal antibody from tissue culture supernatant diluted in phosphate buffered saline, pH 7.4, with protein base, and preserved with sodium azide.

**Applications:**

Anti-Estrogen Receptor antibody reacts with a 67 kDa polypeptide. This antibody stains nuclei of breast epithelial cells and some carcinomas, as well as endometrial epithelia and myometrium. Since the early 1990's the immunohistochemical (IHC) assay determination of Estrogen Receptor status has replaced the dextran-coated charcoal method as a prognostic indicator in breast carcinoma. IHC has shown to be superior in prognostic significance when using any one of several available methods of quantitation using this technique.

**Reactivity:** .....Paraffin, Frozen

**Control:** .....Breast, Endometrium

**Visualization:** .....Nuclear

**Stability:** .....Up to 36 months; store at 2-8° C

**Isotype:** .....IgG<sub>1</sub>/K

Description	Cat. No.	Dilution/ Comments
1 ml, prediluted	CMA660	Ready to use
6 ml, prediluted	CMA661	Ready to use
0.1 ml, concentrate	CMC660-2	1:25-1:100*
0.5 ml, concentrate	CMC661	1:25-1:100*
1 ml, concentrate	CMC662	1:25-1:100*
Positive control slides	CMS603	5 slides per pack

\* The dilutions set forth above are estimates; actual results may differ because of variability in methods and protocols. Validation of antibody performance/protocol is the responsibility of the end user.

*MSDS available upon request.*

**Preparation and Pretreatment:**

1. Cut 3-4 μm section of formalin-fixed paraffin-embedded tissue and place on positively charged slides; dry overnight at 58° C.
2. Deparaffinize, rehydrate, and epitope retrieve; the preferred method is the use of Heat Induced Epitope Retrieval (HIER) techniques using Cell Marque's Trilogy™ in conjunction with a pressure cooker. The preferred method allows for simultaneous deparaffinization, rehydration, and epitope retrieval. Upon completion, rinse with 5 changes of distilled or deionized water.
3. If using HRP detection system, place slides in peroxide block for 10 minutes; rinse. If using AP detection system, omit this step.

**Recommended Protocol for Staining at Room Temperature Using CytoScan™ BSA Detection System:**

1. Apply the antibody and incubate for 60 minutes; rinse.
2. Apply the link and incubate for 10 minutes; rinse.
3. Apply the label and incubate for 10 minutes; rinse.
4. Apply ample amount of chromogen and incubate for 1 to 10 minutes; rinse.
5. Dehydrate and coverslip.

**Recommended Protocol for Staining at Room Temperature Using PolyScan™ Polymer Detection System:**

1. Apply the antibody and incubate for 60 minutes; rinse.
2. Apply the PolyScan™ Polymer Rabbit/Mouse Detection System for 30 min; rinse.
3. Apply ample amount of chromogen and incubate for 1 to 10 minutes; rinse.
4. Dehydrate and coverslip.

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