
# EBV (CS1-4)

 CELL MARQUE D A T A S H E E T

## EBV (CS1-4)

English .....	2	Italiano .....	18
Česky.....	4	Norsk.....	20
Dansk .....	6	Polski.....	22
Nederlands taal .....	8	Português .....	24
Français .....	10	Slovenski .....	26
Deutsch .....	12	Español .....	28
Ελληνικές .....	14	Svensk .....	30
Magyar .....	16		

### References:

1. Murray PG et al. J Pathol. 166: 1-5 (1992)
2. Jarrett RF et al. Blood 78:1-10 (1991)
3. Pailesen G et al. Lancet. 337: 320-322 (1991)
4. Silverberg GS et al Principles and Practice of Surgical Pathology and Cytopathology, 3rd edition. (1997)

**Presentation:**

Anti-EBV is a mouse monoclonal antibody from tissue culture supernatant diluted in phosphate buffered saline, pH 7.4, with protein base, and preserved with sodium azide.

**Applications:**

Anti-Epstein-Barr virus antibody targets the 60 kDa latent membrane protein (LMP-1) encoded by the BNLF1 gene of the Epstein-Barr virus. There is cross-reactivity with Reed Sternberg cells of Hodgkin's disease. The Epstein-Barr virus is an important as a cause of Infectious mononucleosis and has been associated with oral carcinomas.

**Reactivity:** .....Paraffin, Frozen

**Control:** .....Hodgkin's Disease, Infected Tissue

**Visualization:** .....Cytoplasmic

**Stability:** .....Up to 36 months; store at 2-8° C

**Isotype:** .....IgG<sub>1</sub>

Description	Cat. No.	Dilution/ Comments
1 ml, prediluted	CMA424	Ready to use
6 ml, prediluted	CMA425	Ready to use
0.1 ml, concentrate	CMC425-2	1:50 - 1:200*
0.5 ml, concentrate	CMC424	1:50 - 1:200*
1 ml, concentrate	CMC426	1:50 - 1:200*
Positive control slides	CMS424	5 slides per pack

\* The dilutions set forth above are estimates; actual results may differ because of variability in methods and protocols. Validation of antibody performance/protocol is the responsibility of the end user.

*MSDS available upon request.*

**Preparation and Pretreatment:**

1. Cut 3-4 µm section of formalin-fixed paraffin-embedded tissue and place on positively charged slides; dry overnight at 58° C.
2. Deparaffinize, rehydrate, and epitope retrieve; the preferred method is the use of Heat Induced Epitope Retrieval (HIER) techniques using Cell Marque's Trilogy™ in conjunction with a pressure cooker. The preferred method allows for simultaneous deparaffinization, rehydration, and epitope retrieval. Upon completion, rinse with 5 changes of distilled or deionized water.
3. If using HRP detection system, place slides in peroxide block for 10 minutes; rinse. If using AP detection system, omit this step.

**Recommended Protocol for Staining at Room Temperature Using CytoScan™ BSA Detection System:**

1. Apply the antibody and incubate for 60 minutes; rinse.
2. Apply the link and incubate for 10 minutes; rinse.
3. Apply the label and incubate for 10 minutes; rinse.
4. Apply ample amount of chromogen and incubate for 1 to 10 minutes; rinse.
5. Dehydrate and coverslip.

**Recommended Protocol for Staining at Room Temperature Using PolyScan™ Polymer Detection System:**

1. Apply the antibody and incubate for 60 minutes; rinse.
2. Apply the PolyScan™ Polymer Rabbit/Mouse Detection System for 30 min; rinse.
3. Apply ample amount of chromogen and incubate for 1 to 10 minutes; rinse.
4. Dehydrate and coverslip.

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