

Wash Buffers



SAVE UP TO
75%
over leading competitors' brands

Rev. 1.1

Compatible with open, automated IHC instruments.

*TWEEN is a registered trademark of Croda International PLC

Phosphate Buffer Saline (PBS) and **Tris Buffer Saline (TBS)** are concentrated wash solutions designed to rinse reagents from tissue slides. Additionally, they serve to preserve tissue sections between reagent applications. The use of PBS is recommended with any detection system except those that contain alkaline phosphatase. It is recommended that TBS be used when alkaline phosphatase is present to avoid excessive phosphate accumulation and consequent reduction in the rate of reaction. The addition of Tween® 20 to both PBS and TBS may help to reduce background and thereby produce an overall cleaner immunohistochemical stain.

Availability

20x PBS IHC Wash Buffer with Tween® 20	200 ml	934B-06
	1000 ml	934B-09
20x TBS IHC Wash Buffer with Tween® 20	200 ml	935B-06
	1000 ml	935B-09