



unMask hidden pathologies

Melanin Bleach kit

Heavy melanin pigment in the cells can obscure cellular morphology and subsequently impair special stains, IHC, or histochemical staining – and potentially result in an incorrect interpretation. The Melanin Bleach Kit enables the removal of melanin from tissue sections and, in turn, allows for positive staining to be interpreted correctly. It can be used before special stains are applied, or before IHC using DAB.

Benefits

- Save time - no mixing required
- Clean and clear results, easier to analyze
- Reveal cellular morphology obscured by melanin pigment
- Convenient kit format

Learn more at

SigmaAldrich.com/melaninbleachkit

To place an order or receive technical assistance

Order/Customer Service:
SigmaAldrich.com/order

Technical Service:
SigmaAldrich.com/techservice

Safety-related Information:
SigmaAldrich.com/safetycenter

Development/Custom Manufacturing Inquiries:
safcgloba@SigmaAldrich.com

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

Copyright © 2017 EMD Millipore Corporation. All Rights Reserved. MilliporeSigma, and the Vibrant M are trademarks of Merck KGaA, Darmstadt, Germany. Sigma-Aldrich is a trademark of Sigma-Aldrich Co. LLC. or its affiliates. All other trademarks are the property of their respective owners.

Lit. No. RFSHH0021ENUS Ver. 1.0 2017 - 06377 08/2017

**MILLIPORE
SIGMA**